

Rocket

FBS Kits with Sampling Wand

At Rocket Medical we know our customers are our greatest source of new ideas and product improvements. Our latest FBS Kits incorporate your user feedback and will mean better visualisation, easier scalp puncture and sample collection all combined in a smaller tray to save on packaging and storage space.

RELIABLE, COMPATIBLE

ANALYSIS: Rocket FBS Kits are specifically designed for use with Radiometer® ABL™, Siemens® Rapidlab, Instrumentation Laboratory® GEM™, ROCHE OMNI, cobas, Nova Biomedical lactate and other analysers

MAXIMUM ILLUMINATION NOW COMES AS STANDARD:

To give you brighter, more reliable visualisation every kit now contains the **Rocket Amnilume™ Light Source** producing 8800mcd @ 30°.

EASIER SAMPLE COLLECTION:

Losing sight of the puncture wound when transferring between blade and tube holder can be a major problem.

The Rocket Sampling Wand

combines blade and capillary tube holders in one so your view stays clear. It is as easy to use as a ball point pen.



CONTENTS: Endoscope and Rocket Amnilume™ Single Use Light Source, 5 x 10cm x-ray detectable gauze swabs, 6 swab sticks, 2.5gms petroleum jelly, sponge forceps, high heparin capillary tubes, spare capillary tube holder, tube covers, Rocket Sampling Wand
Packed in cartons 10 units, for single patient use.

| Analyser | Minimum recommended sample volume | Rocket® FBS Kit with Sampling Wand | NHS SC Order codes |
|--|-----------------------------------|--|----------------------------|
| ABL™ 5/500, 70, 77 & 90 Series | >65µl | R57024-00-SW | FFY056 |
| ABL™ 800 FLEX with FLEXMODE selected | >35µl | R57012-00-SW R57013-00-SW R57024-00-SW | FFY053 FFY054 FFY056 |
| ABL™ 700 Series with tube volume selected | >35µl | R57018-00-SW | FFY055 |
| Instrumentation Laboratory (IL) GEM™ 4000 Series | 65µl | R57024-00-SW | FFY056 |
| NOVA Biomedical StatStrip® Lactate | 6µl | R57018-00-SW | FFY055 |
| ROCHE® OMNI c, OMNI s cobas b 121, 123, 221 systems | 68-85µl | R57024-00-SW | FFY056 |
| SIEMENS® (Bayer) Rapidlab™ Series | 100µl | R57026-00-SW | FFY057 |

